

Under-Counter Refrigerators UCR 6660 Series



Under-Counter Refrigerators

Under-Counter Refrigerators UCR 6660 Series

Under-Counter Refrigerators UCR 6660 series is equipped with LED digital display and secop compressors. The temperature range of these refrigerators is +2°C~+8°C and comes in a wide range of different no. of shelves and doors. These are easy to operate, easy to maintain and can exhibit better efficiency.

Features

- ◆ R134a is the refrigerant which is free from CFC and foam
- ◆ Electronic thermostat and LED digital display for accuracy and easy reading
- ◆ 100% copper condenser and evaporator
- ◆ 20FT/40HQ are 20ft and 40ft refrigerate containers
- ◆ Auto defrosting and auto evaporation or condensation of water by means of hot pipe in stainless steel box means no floor drain required
- ◆ The internal and external structure is of stainless steel with rounded edges for easy sanitation
- ◆ Insulation 50mm thickness for better energy efficiency
- ◆ Self-closing door by means of spring-loaded hinges ensures fast recovery and energy saving
- ◆ Magnetic-Gasket with high quality neoprene rubber can be replaced without any tools Heating element concealed in door frame to eliminate any condensation
- ◆ Shelf and shelf support are designed to accommodate 1 number of GN 1/1 pan size
- ◆ The doors are structured in offset profile for easy gripping and opening as well as easy maintenance
- ◆ Ventilated refrigeration compartment gives better air circulation, anticorrosion-treated evaporator

Under-Counter Refrigerators

Applications _____

Used mainly in General Lab storage, blood storage, hazardous material storage, Schools & Administrative office, hyper markets and others.

Specification _____

Model No.	UCR 6660	UCR 6661	UCR 6662
Capacity	400L	400L	400L
Temperature Range	+2°C~+8°C	+2°C~+8°C	+2°C~+8°C
No. of Doors	3	3	3
No. of 20FT/40HQ	18/36	18/36	18/36
No. of Drawers	0	2	1
Adjustable Shelves	3	3	3
Product Dimensions	1795*700*850 mm	1795*700*850 mm	1795*700*850 mm
Package Dimensions	1845*740*910 mm	1845*740*910 mm	1845*740*910 mm
Net Weight	115 kg	124 kg	125 kg
Gross Weight	127 kg	140 kg	141 kg
Power Supply	220V/50Hz or 110V/60Hz	220V/50Hz or 110V/60Hz	220V/50Hz or 110V/60Hz

Under-Counter Refrigerators

Model No.	UCR 6663	UCR 6664	UCR 6665
Capacity	400L	400L	400L
Temperature Range	+2°C~+8°C	+2°C~+8°C	+2°C~+8°C
No. of Doors	3	3	3
No. of 20FT/40HQ	18/36	18/36	18/36
No. of Drawers	3	4	6
Adjustable Shelves	3	3	3
Product Dimensions	1795*700*850 mm	1795*700*850 mm	1795*700*850 mm
Package Dimensions	1845*740*910 mm	1845*740*910 mm	1845*740*910 mm
Net Weight	127 kg	140 kg	143 kg
Gross Weight	143 kg	175 kg	159 kg
Power Supply	220V/50Hz or 110V/60Hz	220V/50Hz or 110V/60Hz	220V/50Hz or 110V/60Hz

Optional _____

- ◆ Set drawers instead of doors
- ◆ Castors
- ◆ Stainless steel shelf or plastic-coated shelf is available