

# Under-Counter Refrigerators UCR 7660 Series



# Under-Counter Refrigerators

## Under-Counter Refrigerators UCR 7660 Series

---

Under-Counter Refrigerators UCR 7660 series is equipped with LED digital display and secop compressors. The temperature range of these refrigerators is +2°C~+8°C and comes in a wide range of different no. of shelves and doors. These are easy to operate, easy to maintain and can exhibit better efficiency.

## Features

---

- ◆ R134a is the refrigerant which is free from CFC and foam
- ◆ Electronic thermostat and LED digital display for accuracy and easy reading
- ◆ 100% copper condenser and evaporator
- ◆ 20FT/40HQ are 20ft and 40ft refrigerate containers
- ◆ Auto defrosting and auto evaporation or condensation of water by means of hot pipe in stainless steel box means no floor drain required
- ◆ The internal and external structure is of stainless steel with rounded edges for easy sanitation
- ◆ Insulation 50mm thickness for better energy efficiency
- ◆ Self-closing door by means of spring-loaded hinges ensures fast recovery and energy saving
- ◆ Magnetic-Gasket with high quality neoprene rubber can be replaced without any tools Heating element concealed in door frame to eliminate any condensation
- ◆ Shelf and shelf support are designed to accommodate 1 number of GN 1/1 pan size
- ◆ The doors are structured in offset profile for easy gripping and opening as well as easy maintenance
- ◆ Ventilated refrigeration compartment gives better air circulation, anticorrosion-treated evaporator



# Under-Counter Refrigerators

## Applications

Used mainly in General Lab storage, blood storage, hazardous material storage, Schools & Administrative office, hyper markets and others.

## Specification

Model No.	UCR 7660	UCR 7661	UCR 7662
Capacity	540L	540L	540L
Temperature Range	+2°C~+8°C	+2°C~+8°C	+2°C~+8°C
No. of Doors	4	4	4
No. of 20FT/40HQ	12/30	12/30	12/30
No. of Drawers	0	2	1
Adjustable Shelves	4	4	4
Product Dimensions	2230*700*850 mm	2230*700*850 mm	2230*700*850 mm
Package Dimensions	2280*740*910 mm	2280*740*910 mm	2280*740*910 mm
Net Weight	141 kg	154 kg	155 kg
Gross Weight	161 kg	174 kg	175 kg
Power Supply	220V/50Hz or 110V/60Hz	220V/50Hz or 110V/60Hz	220V/50Hz or 110V/60Hz

# Under-Counter Refrigerators

Model No.	UCR 7663	UCR 7664	UCR 7665
Capacity	540L	540L	540L
Temperature Range	+2°C~+8°C	+2°C~+8°C	+2°C~+8°C
No. of Doors	4	4	4
No. of 20FT/40HQ	12/30	12/30	12/30
No. of Drawers	3	4	6
Adjustable Shelves	4	4	4
Product Dimensions	2230*700*850 mm	2230*700*850 mm	2230*700*850 mm
Package Dimensions	2280*740*910 mm	2280*740*910 mm	2280*740*910 mm
Net Weight	157 kg	167 kg	173 kg
Gross Weight	177 kg	187 kg	193 kg
Power Supply	220V/50Hz or 110V/60Hz	220V/50Hz or 110V/60Hz	220V/50Hz or 110V/60Hz

Optional \_\_\_\_\_

- ◆ Set drawers instead of doors
- ◆ Castors
- ◆ Stainless steel shelf or plastic-coated shelf is available